

CAT URINARY CATHETER

Features:

- One catheter replaces two. Same catheter used to unblock and leave indwelling to drain the bladder. This single device technique is designed to reduce the risk of infection and trauma associated with a two catheter technique.
- FLUSHING STYLET (optional) enables flushing through the stylet during insertion.
- Moveable new suture device allows length adjustment.
- Non-irritating radiopaque polyurethane material.
- Rigid for easy insertion.
- Softens indwelling at body temperature.
- Finished open tip.

Sizes & Packings

- FG 3 (OD: 1.00 mm & Length: 110 mm)
- FG 4 (OD: 1.30 mm & Length: 110 mm)

